




## FITTING INSTRUCTIONS FOR THE KNEE EXTENSION

1. Place the patient in a reclined or resting position.
2. ONLY open BACK CUFFS #1 and #2.
3. Place directly over KNEE with 2" strap above and 2" strap below KNEE CAP.
4. Open and tighten straps above and below KNEE CAP until cams align with center of knee joint.
5. Close cuff #1 (BACK OF THIGH) and pull to line.
6. Bend splint at hinge pushing toward the leg.
7. Close cuff #2 (CALF) and pull to line.
8. Support the leg with pillow to elevate the foot and float the heel.

OPTIONAL: If MPO Boot is being worn in conjunction with the Knee Extension Dynasplint® System, apply the boot first following the enclosed instructions in packaging.



**DO NOT OPEN CUFFS THAT HAVE BEEN SECURED WITH DUCT TAPE OR GLUE**